

1.本地运行neuroglancer


(1) 从github上下载代码：

安装git (<https://baijiahao.baidu.com/s?id=1601036689157983619&wfr=spider&for=pc>) ,

创建GitHub账号，配置GitHub和本机电脑的联系 (https://blog.csdn.net/weixin_43788566/article/details/84892542?ops_request_misc=%257B%2522request%255Fid%2522%253A%2522163394855316780265439001%2522%252C%2522scm%2522%253A%25220140713.130102334.pc%255Fall.%2522%257D&request_id=163394855316780265439001&biz_id=0&utm_medium=distribute.pc_search_result.none-task-blog-2~all~first_rank_ecpm_v1~rank_v31_ecpm-1-84892542.first_rank_v2_pc_rank_v29&utm_term=github%E7%9A%84%E6%9C%AC%E5%9C%B0%E9%85%8D%E7%BD%AE&spm=1018.2226.3001.4187) 。

找到neuroglancer， (<https://github.com/google/neuroglancer>)

找一文件夹，在命令行输入 `git clone git@github.com:google/neuroglancer.git`



MINGW32:/d/WorkSpace

```
stephen@LAPTOP-3PQFLQQN MINGW32 /d/WorkSpace
$ git clone git@github.com:google/neuroglancer.git
Cloning into 'neuroglancer'...
remote: Enumerating objects: 13591, done.
remote: Counting objects: 100% (109/109), done.
remote: Compressing objects: 100% (78/78), done.
remote: Total 13591 (delta 56), reused 48 (delta 31), pack-reused 13482
Receiving objects: 100% (13591/13591), 6.31 MiB | 2.08 MiB/s, done.
Resolving deltas: 100% (10427/10427), done.
```

```
stephen@LAPTOP-3PQFLQQN MINGW32 /d/WorkSpace
$ |
```

(2) 本地运行neuroglancer

找到readme，按照步骤执行，执行完第四步，本地就可以运行起来了，建议下载个vscode，在其命令行上执行以下指令。

Building

node.js is required to build the viewer.

1. First install NVM (node version manager) per the instructions here:

<https://github.com/creationix/nvm>

2. Install a recent version of Node.js if you haven't already done so:

```
nvm install stable
```

3. Install the dependencies required by this project:

(From within this directory)

```
npm i
```

Also re-run this any time the dependencies listed in [package.json](#) may have changed, such as after checking out a different revision or pulling changes.

4. To run a local server for development purposes:

```
npm run dev-server
```

This will start a server on <http://localhost:8080>.

5. To run the unit test suite on Chrome:

```
npm test
```

To run only tests in files matching a given glob pattern:


```
npm test -- --pattern='<pattern>'
```

For example,

```
npm test -- --pattern='src/neuroglancer/util/uint64*'
```

6. See [package.json](#) for other commands available.

最后一条运行指令执行完，在本地浏览器 输入 <http://127.0.0.1:61849>






2.搭建本地文件服务器









经过多天调研，我选用的是http-server这款web服务器框架，

安装 (https://blog.csdn.net/weixin_37856861/article/details/83274047?ops_request_misc=%257B%2522request%255Fid%2522%253A%2522163394966716780265411404%2522%252C%2522scm%2522%253A%25220140713.130102334.pc%255Fall.%2522%257D&request_id=163394966716780265411404&biz_id=0&utm_medium=distribute.pc_search_result.none-task-blog-2~all~first_rank_ecpm_v1~rank_v31_ecpm-6-83274047.first_rank_v2_pc_rank_v29&utm_term=windows+%E5%AE%89%E8%A3%85+http-server&spm=1018.2226.3001.4187)

安装成功后，文件服务器搭建成功，



Index of /

 (drw-rw-rw-)	1280um/
 (drw-rw-rw-)	160um/
 (drw-rw-rw-)	20um/
 (drw-rw-rw-)	320um/
 (drw-rw-rw-)	40um/
 (drw-rw-rw-)	640um/
 (drw-rw-rw-)	80um/
 (-rw-rw-rw-) 1.2k	info

Node.js v14.17.6/ [http-server](#) server running @ 10.100.159.24:8080

3.转化数据为neuroglancer所需格式

(1) 安装“转换脚本项目”的环境

我在GitHub上找到一转化格式的脚本 (<https://github.com/HumanBrainProject/neuroglancer-scripts>)，修改后既能生成neuroglancer所需数据源。

按照步骤安装此项目，运行此脚本前，需进入其自带的虚拟环境。note：红框处命令必须执行

neuroglancer-scripts

Tools for converting volumetric images and surface meshes to the pre-computed format of [Neuroglancer](#).

pypi
v0.3.0
build
passing
codecov
92%
docs
passing

Installation

The easiest way to install the latest stable version of neuroglancer-scripts is through `pip`. Using a virtual environment is recommended:

```
python3 -m venv venv/
. venv/bin/activate
pip install neuroglancer-scripts
```



Usage

See the [documentation](#).

Development

The code is hosted on <https://github.com/HumanBrainProject/neuroglancer-scripts>.

Useful commands for development:

```
git clone https://github.com/HumanBrainProject/neuroglancer-scripts.git

# Install in a virtual environment
cd neuroglancer-scripts
python3 -m venv venv/
. venv/bin/activate
pip install -e .[dev]

# Tests
pytest # run tests
pytest --cov=neuroglancer_scripts --cov-report=html # detailed test coverage report
```

(2) 根据步骤转换格式

进入这个网站 (<https://neuroglancer-scripts.readthedocs.io/en/latest/examples.html>)

执行以下4步 (也可按照下面我详解的步骤操作) :

Conversion of BigBrain

BigBrain is a very large image ($6572 \times 7404 \times 5711$ voxels) reconstructed from 7404 serial coronal section of a human brain, with a resolution of about 20 microns.

1. Download slices from

ftp://bigbrain.loris.ca/BigBrainRelease.2015/2D_Final_Sections/Coronal/Png/Full_R

2. Create `info_fullres.json` with the appropriate metadata:

```
{
  "type": "image",
  "data_type": "uint8",
  "num_channels": 1,
  "scales": [
    {
      "chunk_sizes": [],
      "encoding": "raw",
      "key": "full",
      "resolution": [21166.666666666666, 20000, 21166.666666666666],
      "size": [6572, 7404, 5711],
      "voxel_offset": [0, 0, 0]
    }
  ]
}
```

3. Create raw chunks

```
generate-scales-info info_fullres.json 8bit/
slices-to-precomputed --input-orientation RIA <path/to/slices> 8bit/
compute-scales --outside-value=255 8bit/
```

4. Optionally, convert raw chunks to JPEG:

```
generate-scales-info --encoding=jpeg 8bit/info jpeg/
convert-chunks --jpeg-plane=xz 8bit/ jpeg/
```

详解步骤:

重新拉取我github上的neuroglancer-scripts项目 (<https://github.com/stephen-zhouyang/neuroglancer-scripts>)，我对其做了修改，建议用pycharm打开，：

①安装同上：

按照步骤安装此项目，运行此脚本前，需进入其自带的虚拟环境。note：**红框处命令必须执行**

neuroglancer-scripts

Tools for converting volumetric images and surface meshes to the pre-computed format of [Neuroglancer](#).

pypi
v0.3.0
build
passing
codecov
92%
docs
passing

Installation

The easiest way to install the latest stable version of neuroglancer-scripts is through `pip`. Using a virtual environment is recommended:

```
python3 -m venv venv/
. venv/bin/activate
pip install neuroglancer-scripts
```



Usage

See the [documentation](#).

Development

The code is hosted on <https://github.com/HumanBrainProject/neuroglancer-scripts>.

Useful commands for development:

```
git clone https://github.com/HumanBrainProject/neuroglancer-scripts.git

# Install in a virtual environment
cd neuroglancer-scripts
python3 -m venv venv/
. venv/bin/activate
pip install -e .[dev]

# Tests
pytest # run tests
pytest --cov=neuroglancer_scripts --cov-report=html # detailed test coverage report
```

```
. venv/bin/activate
```


这条指令的含义是进入其自带虚拟环境，我们也可以这样做，当导入pycharm后，后自动生成venv library root，然后进入目录，执行

```
.\activate
```

即进入其虚拟环境，再执行

```
pip install neuroglancer-scripts
```

虚拟环境如下：




```
nib-conform.exe
Terminal: Local x +
origin git@github.com:stephen-zhouyang/neuroglancer-scripts.git (push)
(venv) D:\WorkSpace\otherWorkSpace\neuroglancer-scripts\src\neuroglancer_scripts\scripts>
(venv) D:\WorkSpace\otherWorkSpace\neuroglancer-scripts\src\neuroglancer_scripts\scripts>
(venv) D:\WorkSpace\otherWorkSpace\neuroglancer-scripts\src\neuroglancer_scripts\scripts>
```

②通过ftp工具filezilla (<https://filezilla-project.org/download.php>) 下载源文件:

下载安装好后, 红框输入主机名称, 进入到以下目录:

```
ftp://bigbrain.loris.ca/BigBrainRelease.2015/2D_Final_Sections/Coronal/Png/Full_Resolution/
```



下载名称为pm0001-pm0020的20张图片d到本地。

③然后依次执行以下脚本, 记住更换20张图片所在目录, 我的目录是 D:/datatest/precomputed/png20 8bit/, 更换成你们自己的图片所在的目录。当执行python change8bitName.py和python changeJpegName.py时, 也要更换其中的目录(见下图)。

生成raw型文件

```
generate-scales-info info_fullres.json 8bit/
slices-to-precomputed --input-orientation RIA D:/datatest/precomputed/png20
8bit/
compute-scales --outside-value=255 8bit/
python change8bitName.py
```

生成jpeg文件


```
generate-scales-info --encoding=jpeg 8bit/info jpeg/  
convert-chunks --jpeg-plane=xz 8bit/ jpeg/  
python changeJpegName.py
```

```
changeJpegName.py ×  
1 import os  
2  
3  
4 root = 'D:\WorkSpace\otherWorkSpace\neuroglancer-scripts\src\neuroglancer_scripts\scripts\jpeg'  
5 print(os.listdir(root))  
6  
7 for folder in os.listdir(root):  
8     f = root + "\\" + folder  
9     for s1 in os.listdir(f):  
10         f1 = f + "\\" + s1  
11         for s2 in os.listdir(f1):  
12             f2 = f1 + "\\" + s2  
13             for s3 in os.listdir(f2):  
14                 f3 = f2 + "\\" + s3  
15                 ss = s1 + "_" + s2 + "_" + s3  
16                 open(f + "\\" + ss, 'wb').write(open(f3, 'rb').read())  
17
```

4.调试

1.启动http-server

在生成的jpeg目录下打开命令行，输入以下指令：

```
http-server --cors 127.0.0.1 -p 8080
```

note: --cors是解决跨域问题的。

« neuroglancer-scripts » src » neuroglancer_scripts » scripts » jpeg »			
	名称	修改日期	类型
ds	20um	2021/10/11 21:59	文件夹
	40um	2021/10/11 21:59	文件夹
	80um	2021/10/11 21:59	文件夹
	160um	2021/10/11 21:59	文件夹
	320um	2021/10/11 21:59	文件夹
	640um	2021/10/11 21:59	文件夹
	1280um	2021/10/11 21:59	文件夹
	info	2021/10/11 21:59	文件

成功结果如下，

MINGW32:/d/WorkSpace/otherWorkSpace/neuroglancer-scripts/src/neuroglancer_scripts/scripts/jpeg

```
stephen@LAPTOP-3PQFLQQN MINGW32 /d/WorkSpace/otherWorkSpace/neuroglancer-scripts
/src/neuroglancer_scripts/scripts/jpeg (master)
$ http-server --cors 127.0.0.1 -p 8080
Starting up http-server, serving ./

http-server settings:
CORS: 127.0.0.1
Cache: 3600 seconds
Connection Timeout: 120 seconds
Directory Listings: visible
AutoIndex: visible
Serve GZIP Files: false
Serve Brotli Files: false
Default File Extension: none

Available on:
http://10.100.159.24:8080
http://127.0.0.1:8080
Hit CTRL-C to stop the server
```

在本地启动的neuroglance中输入 `precomputed://http://127.0.0.1:8080` 即完全成功

